

CYP2B6 Polyclonal Antibody



Catalog Number: 162186

Product name

CYP2B6 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CYP2B6

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 160-340 of human CYP2B6 (NP_000758.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000


IHC 1:50 - 1:100

Validations



Western blot - CYP2B6 Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using CYP2B6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - CYP2B6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung using CYP2B6 antibody at dilution of



1:100 (40x lens).


CYP2B6 Polyclonal Antibody



Catalog Number: 162186

Immunohistochemistry - CYP2B6 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human colon using CYP2B6 antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - CYP2B6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using CYP2B6 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - CYP2B6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using CYP2B6 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - CYP2B6 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human stomach using CYP2B6 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - CYP2B6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using CYP2B6 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - CYP2B6 Polyclonal Antibody

CYP2B6 Polyclonal Antibody



Catalog Number: 162186


Immunohistochemistry of paraffin-embedded mouse lung using CYP2B6 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CYP2B6 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse liver using CYP2B6 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CYP2B6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using CYP2B6 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CYP2B6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using CYP2B6 antibody at dilution of 1:100 (40x lens).