

CRYGC antibody



Catalog Number: 109586

Product name

CRYGC antibody

Immunogen

[Human CRYGC Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

CRYGC Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 21kd

Preparation

This antibody was obtained by immunization of CRYGC recombinant protein (Accession Number: NM_020989). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109586(CRYGC antibody) at dilution of 1:500