

ELOF1 antibody



Catalog Number: 110286

Product name

ELOF1 antibody

Immunogen

[Human ELOF1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ELOF1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse pancreas tissue. Observed molecular weight by Western-blot: 9 kDa

Preparation

This antibody was obtained by immunization of ELOF1 recombinant protein (Accession Number: XM_017027359). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110286(ELOF1 antibody) at dilution of 1:500