

# CSNK1A1 Polyclonal Antibody



Catalog Number: 163096

## Product name

CSNK1A1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to CSNK1A1

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CSNK1A1 (NP\_001020276.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - CSNK1A1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CSNK1A1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.