

# Anti-SENP1 Rabbit antibody



Catalog Number: 166896

## Product name

Anti-SENP1 Rabbit antibody

## Specificity

Human

## Antibody description

SENP1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human SENP1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/1000-1/5000

IHC: 1/20-1/100


ICC/IF: 1/20-1/100

FC: 1/20

## Validations



Western blot detection of SENP1 in 293, Hela, COS7, PC3 cell lysates using SENP1 Rabbit pAb(1:1000 diluted). Predicted band size: 74 kDa. Observed band size: 76 kDa.



Immunocytochemistry of SENP1 (green) in Hela cells using SENP1 Rabbit pAb at dilution 1/50, and DAPI (blue)