

YTHDF3-specific antibody



Catalog Number: 116896

Product name

YTHDF3-specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

YTHDF3-specific Rabbit Polyclonal antibody. Positive IF detected in A549 cells. Positive WB detected in HepG2 cells, A549 cells, HeLa cells, HL-60 cells, mouse liver tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of YTHDF3-specific recombinant protein (Accession Number: NM_152758). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116896(YTHDF3 Antibody) at dilution of 1:600



Immunofluorescent analysis of A549 cells using Catalog No:116896(YTHDF3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)