

IL-28A antibody



Catalog Number: 107396

Product name

IL-28A antibody

Specificity

Human; other species not tested.

Antibody description

IL-28A Mouse Monoclonal antibody. Positive WB detected in ag18525. Positive IHC detected in human tonsillitis tissue, human small intestine tissue. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of IL-28A recombinant protein (Accession Number: NM_172138). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:5000-1:50000

IHC: 1:20-1:200


Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107396(IL28A Antibody) at dilution of 1:5000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107396(IL28A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107396(IL28A Antibody) at dilution of 1:50

IL-28A antibody



Catalog Number: 107396

(under 40x lens)